






**Accessorie:**  
**Mounting kit / Art.-no.: 164.30.00**  
 Walking stick holder for rollator: Polo-Plus-T, Jumbo






**REBOTEC®** Rehabilitationsmittel GmbH  
 DE-49610 Quakenbrück, Artlandstr. 57-59  
 Telephone: +49(0) 5431/9416-0  
 Telefax : +49(0)5431/9416-66  
<http://www.rebotec.de>  
 e-mail: info@rebotec.de

MA-Nr.: 02.10.140/2014-02

<p><b>Mounting kit</b></p> <p><b>A. clamp holder</b> D=25 mm</p>  <p><b>B. quiver holder</b> D=30mm</p> 	<p><b>2 versions: mounting and adjustment</b></p> 
---	---

<p><b>Mounting (same assembly operation of the parts A and B</b></p> 	<p><b>clamp and screw the holder</b></p> 
--	--

<p><b>Position rotatable 360 degrees</b></p> 	<p><b>Situating the walking stick</b></p> 	<p><b>Clamping the walking stick</b></p> 
--	---	--